

Product Name :
Pocket Concrete Penetrometer

Product Code :
NLAB -CONCRLAB120034



Description :

Pocket Concrete Penetrometer

Technical Specification :

Scale range is 0-50 kg/cm² when the penetration resistance reaches a value of 35kg/cm² the concrete is assumed initially set. Supplied complete in carrying case.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003