

Product Name :
Pocket Microscope

Product Code :
NLAB-EDUSCHOOLEQ173077

Description :

Pocket Microscope

Technical Specification :

Portable LED lighted Pocket Microscope.

Runs on 3 SG3 Button Cell Batteries (Included).

The microscope also includes a built-in LED light that provides a bright, clear image, along with a rubberized eyepiece for comfortable viewing.

Measures 1.25 x 3.25 x 0.75 inches (W x H x D).

Power and zoom controls (60x-75x magnification) for bright, clear images.

So lightweight you can bring it anywhere.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003