

**Product Name :**  
Pocket Penetrometers

**Product Code :**  
NLAB-SOILELAB140023



**Description :**

Pocket Penetrometers

**Technical Specification :**

Designed for making field classification of cohesive soils in terms of consistency, shear strength and approximate unconfined compressive strength.

Stainless steel construction.

Tip Diameters: 4.5 mm dia. for very hard soil, 6.35 mm for medium and soft soil, 8.98 mm for soft soil.

Measuring range: 0 to 1000 kPa

Dimensions (assembled): 210 mm length x 20 mm dia. approx.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

