

**Product Name :**  
Portable Cylinder Boring Bar

**Product Code :**  
NLAB-MACHINELABWK32009



**Description :**

Portable Cylinder Boring Bar

**Technical Specification :**

**SPECIAL FEATURES :**

- High Precision
- No Need of Honing
- High Efficiency
- Adopton of Guiders
- Setting by Special Micrometer
- Easily Portable

**Standard Equipment :**

Single Phase Motor A/c 220 / 250V with flexible wire & 3 pin

plug, Equipment Box, Micrometer, Tool Lapping Fixture, Lapping Disc, Tool Holders with one cutter bit, Anchoe Assemble & 5 sets of Guiders (for 001 3 set of Guiders) & Necessary Operating Spanners.

~~0040~~" to 2.200"

~~0090~~" to 3.000"

~~0090~~" to 4.200"

~~0070~~ Single Speed

~~0050~~ Double Speed

~~0066~~ to 7.000"

~~0080~~" to 9.600"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003