

Product Name :
Portable Digital Hardness Tester

Product Code :
TN185



Description :

Portable Digital Hardness Tester

Technical Specification :

Portable Digital Hardness Tester

Features:

Dynamic test indicator

Large digital readout

Hardness values in all major international scales with simple conversion facilities from one to the other

RS232 output for connection to PC or serial printer

Memory storage range for in excess of 400 readings

Statistical summary

Integral icon facility provides operator with easy visual identification of mode in which unit is operating

Operates in temperatures from -5°C to +35°C

Upper and lower control limits

Last reading recall

Supplied with two certified test blocks, adjuster key and carrying case

Technical Specification:

Hardness Scale Vickers Pyramid No., Brinell, Rockwell B, Rockwell C,

Tensile Strength, Shore Scleroscope

Testing Range Vickers Pyramid No. 35-1000

Brinell 100-500

Rockwell B 30-100

Rockwell C 20-70
Tensile Strength 255-1999 N/mm²
Shore Scleroscope 24-97
Resolution 1 Vickers Pyramid No., Brinell, Tensile Strength, Shore Scleroscope
0.1 Rockwell B, Rockwell
Power 9V battery or adaptor
Output RS-232 serial output

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003