

Product Name :
Portable Electronic Thermo-Hygrometer

Product Code :
NLAB-EDUCATIONTVET15030



Description :

Portable Electronic Thermo-Hygrometer

Technical Specification :

Combined HR probe with integrated temperature sensor, cable 1 m

Ranges:

Temperature: 0.0 to 60.0 ° C

RH: 20.0 to 95.0%

Precision :

Temperature: +/- 0.5 ° C

RH: 3% RH (50 to 85% RH), +/- 4% RH (outside)

Power supply: 1 x 9 V (auto-extinction)

Supplied in a box with a probe and battery

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003