

Product Name :
Portable Steam Autoclave (Semi Automatic)

Product Code :
NLAB-LABORATORYEM47008



Description :

Portable Steam Autoclave (Semi Automatic)

Technical Specification :

Portable Steam Autoclave (Semi Automatic) -

Features:-

Fully Stainless Steel Construction

Stainless Steel Seamless diepressed deep drawn chamber.

Unique Shape heater of outer SS tubing.

Over heat Auto shut down at 1370 C.

?p based PID controller with dual display for Temp. & Time

R & D Centers

Stainless Steel Safety Valve

Microprocessor Digital temperature controller cum timer

Dual Display

Specially designed heater for optimum chamber usage

Color coded monitoring gauge

Vacuum Breaker to remove air packets in the chamber

Food grade silicon rubber gasket for better durability

Calibrated temperature cum pressure gauge

Specifications:-



Capacity
0.1 lbs

Dimension
200 x 200mm

Power Consumption
2.5 Kw

0.22 atkg/cm2
pressure

0.10 C/1340 C
Temperature

6-10 min
time

230V
Voltage

65 kg
Weight

We are leading manufacturers, suppliers of Portable Steam Autoclave (Semi Automatic) for Autoclave. Contact us to get high quality Portable Steam Autoclave (Semi Automatic) for Autoclave for schools, colleges, universities, research labs, laboratories and various industries.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003