

Product Name :
Portable Tissue Culture Rack

Product Code :
NLABMEDICAL-M3160002

Description :

Portable Tissue Culture Rack

Technical Specification :

Tissue culture rack are very useful for various types of tissue culture plants and other tissue culture work in tissue culture room.

Made By: Tissue rack are made of tubular mild steel pipes With Epoxy Powder Coated.

Racks Specification : Each rack of size 48"x21" (depth) is covered with thick glass or Acrylic. Distance between each rack is 15".

Lighting : Four fluorescent tubes are glare free lights are built in upper portion in worktable.

UV Germicidal Light : Are fitted inside the rack.

Chocks : All fluorescent tubes chocks and ultra violet germicidal tube chocks are fitted out side of rack.

Timer : 0-24 Hours Digital Timers are provided of tissue culture rack.

Control Panels : All control panels are fitted tops of tissue culture racks. On/Off switch of fluorescent tubes and UV tubes.

All control panels are fitted separately in each rack with indicator light.

Easy Mobility: The whole rack is provided with four nos. castor wheels at the bottom for easy mobility to any place.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

