

**Product Name :**  
Positioning Electronic Total Station

**Product Code :**  
NLAB-ENGINEERINGLB85010



**Description :**

Positioning Electronic Total Station

**Technical Specification :**

Battery

Battery charger

Plastic rain cover and instruction manual

Magnification: 30 X

Least Count : 1"

Distance Measuring Range With Single

Accuracy: 2"

Prism: 2000 M

Distance Accuracy : +(2mm + 2 ppm) XD

Toolkit with case

Plastic carrying case

Field of view: 1o30'

Display: Two Side (Graphics LCD - 24 Alpha Num Keys)

Internal Memory: 24000 Points

Downloading software with F-4 data transfer cable 1 No.

Wooden Tripod 1 No

Prism Holder SIN-Tilt 3 1 No.

Single Prism Assembly: Single Prism 1 No.

Prism Pole with bubble 1 No.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003