

Product Name :
Pounding Tester Machines

Product Code :
NLAB-TESTINGLABM53001



Description :

Pounding Tester Machines

Technical Specification :

Pounding Tester - We offer Pounding Tester, which is used to test the compression of foam. The effect of repeated application of pressure on the surface of a slab of polyurethane foam experienced in its normal life is evaluated by simulating these conditions in laboratory by applying and releasing a specified force for a specified number of times, after which the effect of such pounding on the physical characteristics such as hardness index, thickness, and cell structure are observed to evaluate the quality of foam slab.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003