

Product Name :
Power LED Head Light

Product Code :
NLABMEDICAL-M3480002

Description :

Power LED Head Light India, Power LED Head Light manufacturer India, suppliers, Medical Power LED Head Light India, Hospital Power LED Head Light India

Technical Specification :

Light Source: Super-bright and low light attenuation.

Dimensions: Head 1.06 inches.

Housing Color: Black.

Switch: Button switch on the light body of headlamp (Three different modes).

Three Modes: High beam, low beam and strobe light.

Focusing and Zoom: Stretching out for focusing and drawing back for zoom.

Lens: Convex optic lens, good spot lighting and uniform light.

Battery Powered: A piece of 18650 lithium ion battery (Included).

Charging Method: Direct battery charging.

Continuous Run Time: 2-3 Hours (Based on different modes).

Light Color: Neutral white.

Luminous Flux: Up to 280-350 lm at different modes.

Illumination Distance: Within 600 feet on high beam.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments