

Product Name :
Power Supply

Product Code :
NLAB-EDUSCHOOLEQ45040



Description :

Power Supply

Technical Specification :

This power source can be used on 75V or 150V with a maximum current of 250 milliamps.

Run time for gels: 40-50 minutes at 75V; 20-30 minutes at 150V.

The power supply can be used to run one cell or two cells simultaneously.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003