

**Product Name :**  
Power Transistor Characteristics Trainer Digital

**Product Code :**  
NLAB-ELECTRICALAB210041



**Description :**

Power Transistor Characteristics Trainer Digital

**Technical Specification :**

Digital Meters:

Voltmeter 20V DC.

Ammeter 200mA DC.

Voltmeter 200V DC.

Power supplies:

DC Supply IC Regulated 0-5V DC, 150mA.

DC Supply IC Regulated 0-30V DC, 150mA.

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Components are mounted on the panels are:

Power Transistor 2N3055

Voltage Control through Potentiometer.

Salient Features:

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003