

Product Name :
PPR Welding Machine

Product Code :
NLAB-VOCATIONALTVET16011



Description :

PPR Welding Machine

Technical Specification :

Size:-

20mm to 63mm, Copper head.

Temperature regulator to adjust the welding temp, Two heating switches.

Includes:-

Pipe cutter.

Measuring tape.

Allen key.

Holding stand.

Carrying metal case

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003