

**Product Name :**  
Practical LAN Network Trainer

**Product Code :**  
NLAB-ENGINEERINGLB54055



**Description :**

Practical LAN Network Trainer

**Technical Specification :**

This trainer provides essential practical experiments for overall program in meeting A+ Certification and Network Technology Practices as well as Network+ requirements.

The LAN Trainer is developed to give students practical as well as theoretical skills in setting up network systems from simple to complex.

It is a modular system and can start at a low level and work up to more complex LAN & WAN systems.

We are well-known manufacturers, OEM suppliers of Practical LAN Network Trainer for Applied Machines. Contact us for high quality Practical LAN Network Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003