


Product Name :
Pre Pinch Cum Pyrmid Type Plate Rolling Machine

Product Code :
NLAB-MACHINELABWK16007



Description :

Pre Pinch Cum Pyrmid Type Plate Rolling Machine

Technical Specification :

Pre Pinch Cum Pyrmid Type Plate Rolling Machine .

- Three Roll Symmetrical Principle for Pre-pinching of Plate ends & Rolling Without Removing the Plate.
- Modern Design, Robust Steel Construction.
- High Quality Alloy Steel Rolls.
- Built in Cone Rolling Feature.
- Precise Controlls Through Push Button Operation with Interlock Circuit.
- Electrical / Mechanical Overload Protection.
- Centralised Lubrication System.

- Rool Position Setting by Dial / Scale Readout.

CAPACITY RANGE :

- Plate Thickness : 5mm to 50mm
- Plate Width : 1000mm to 3000mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003