


Product Name :
Precision Balance, Type Pa 4102 C

Product Code :
NLAB-EDUCATIONTVET12010



Description :

Precision Balance, Type Pa 4102 C

Technical Specification :

Application: weighing (18 + 1 units), counting, percentage weighing

LCD display screen

Communication interface: RS-232

Spirit level, leveling feet

Weighing hook under scale

Stainless steel weighing pan: Ø 180 mm

Internal calibration

Weight: 4.100 g

Resolution: 0.01 g

Power supply: 220V - 50Hz

Weighing hook under scale

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003