


Product Name :
Precision Balance Type PR 4202 / E

Product Code :
NLAB-EDUCATIONTVET13011



Description :

Precision Balance Type PR 4202 / E

Technical Specification :

Basic weighing (10 units), piece counting, weighing in%

Backlit liquid crystal display (LCD)

Communication interface: RS-232

Spirit level, leveling feet

Weighing hook under scale

Stainless steel weighing pan: Ø 180 mm

External calibration

Weight: 4.200 g

Display resolution: 0.01 g

Repeatability: 0.01 g

Power supply: 230V - 50Hz

Calibration mass, class F1, 2 kg, delivered with certificate

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003