

Product Name :
Precision Balance

Product Code :
NLAB-EDUCATIONTVET12009



Description :

Precision Balance, Type Pa 413 C

Technical Specification :

Application: weighing (18 + 1 units), counting, percentage weighing

LCD display screen

Communication interface: RS-232

Spirit level, leveling feet

Weighing chamber with 3 transparent sliding doors

Weighing hook under scale

Stainless steel weighing pan: Ø 120 mm

Internal calibration

Weight: 410 g

Resolution: 1 mg

Power supply: 220V - 50Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003