

Product Name :
Prescription Balance

Product Code :
NLAB-SCIENTIFICEQ28017



Description :

Prescription Balance

Technical Specification :

Preferably Torsion with a removable glass pan or S.S pan suitable for weighing 50gms to 200gms and sensitivity 30gms for one division indicator shifting on the index plates (Known as dispensing balance)

Available in following ranges:

50gms

100gms

200gms

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

