

Product Name :
Pressure Flow Packed Bed

Product Code :
NLAB-TECHNICALAB32005



Description :

Pressure Flow Packed Bed

Technical Specification :

We are one of the leading manufacturers, suppliers and exporters of optimum quality Pressure Flow Packed Bed. This equipment is highly appreciated among our customer for their flawless design, hassle free installation and easy maintenance. The equipment offered by us is equipped with a rotameter and a manometer for the flow and drop measurement of water respectively.

Electricity 500 watts, 220 V, 1 Phase

Floor Area 0.75 m x 0.75 m

Packed Column : Material Borosilicate Glass with both end made of Stainless Steel

Packing : Material Borosilicate Glass Raschig rings, Size 8-10 mm approx

Water tank : Material Stainless Steel, Capacity 30 Litres

Water Circulation : FHP Pump Crompton/Sharp make

Flow Measurement : Rotameter for water

Drop Measurement : Manometer

The whole set-up is well designed and arranged in a good quality painted structure

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com