

**Product Name :**  
Pressure Measurement Trainer

**Product Code :**  
NLAB-ELECTRONICSAB110005



**Description :**

Pressure Measurement Trainer

**Technical Specification :**

Pressure Measurement Trainer Lab Equipments Manufacturer, Suppliers & Exporters PRESSURE MEASUREMENT TRAINER

To study of Pressure

Range :- 300mm.

Pressure Cell

Sockets at different places for observing / measuring the signals

3.5 digit DPM indication of Pressure

220V, Built in IC regulated power supply

Enclosed in an attractive ABS plastic cabinet with cover for safety of components.

Pressure generator with dial display.

Detailed instruction manual

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

