

Product Name :
Pressure Module

Product Code :
NLAB-LABORATORYEM26010



Description :

Pressure Module

Technical Specification :

AC integrated class 100 integration pressurizing module to pressurize and sustain the clean air with HEPA filter of 0.3 micron, pre filter, motor blower assembly. It delivers cooled and filtered air into the room.

In return air duct having one pre filter of 5 micron, 5 % fresh air provision with dampers for control the fresh air quantity.

Technical Specification:

Air changes / Hour

40 Temp

24C Humidity

45-60% Class

10,000 Noise level

65 decibels +/-5 Standard

FED STD 209E

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003