

Product Name :
Pressure Transducer Trainer

Product Code :
NLAB-ELECTRONICSAB110035



Description :

Pressure Transducer Trainer

Technical Specification :

SCOPE OF LEARNING:

Study of Pressure Measurement using Pressure Cell

TECHNICAL SPECIFICATIONS:

Digital Meters:

Digital Voltmeter : 0 to 10V

Power Supplies:

DC Supply : Built in IC regulated power supplies

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Range :- 0 -10 KG/300mm.

Pressure Cell

Pressure generator with dial display

Pressure Valve -01 no. (manually)

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x10"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003