

<b>Product Name :</b> Prism Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173052
<b>Description :</b>	
Prism Microscope	
<b>Technical Specification :</b>	
Eyepieces: Wide Field 10x	
Total Magnification: 40x ~ 400x	
Conjugated Distance: 145 mm.	
Stand: Heavy base stand made of die cast metal for the body stability.	
Focusing: By moving stage, Range : 10 mm.	
Body: Stable and robust, assembled from aluminium metal parts, Monocular tube inclinable upto 90° mechanical tube length 130 mm. Overall height of the microscope is 285 mm.	
Illumination: Through prism.	
Nosepiece: triple revolving nosepiece with positive centering and click stops.	
Accessories: Vinyl Cover, Duster, and Operating Manual.	
Achromatic Objectives: 4x, 10x and 40x	
Working Stage: Fixed square metal stage Size 100 x 100mm. with two slide holding clips.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003