

**Product Name :**  
Prism Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ173052

**Description :**

Prism Microscope

**Technical Specification :**

Eyepieces: Wide Field 10x

Total Magnification: 40x ~ 400x

Conjugated Distance: 145 mm.

Stand: Heavy base stand made of die cast metal for the body stability.

Focusing: By moving stage, Range : 10 mm.

Body: Stable and robust, assembled from aluminium metal parts, Monocular tube inclinable upto 90° mechanical tube length 130 mm. Overall height of the microscope is 285 mm.

Illumination: Through prism.

Nosepiece: triple revolving nosepiece with positive centering and click stops.

Accessories: Vinyl Cover, Duster, and Operating Manual.

Achromatic Objectives: 4x, 10x and 40x

Working Stage: Fixed square metal stage Size 100 x 100mm. with two slide holding clips.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003