

Product Name :
Prism Mould 75x75x285 mm to UNI

Product Code :
NLAB-CEMETLAB130023



Description :

Prism Mould 75x75x285 mm to UNI

Technical Specification :

Two gang prism mould 75x75x285 mm complete with contact points.

Gauge length 250mm, total length 285 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003