

**Product Name :**  
Prismatic Compass

**Product Code :**  
NLAB-ENGINEERINGLB90001



**Description :**

Prismatic Compass

**Technical Specification :**

This precision engineered Prismatic Compass is one of the important devices which imparts immensely accurate readings for navigational purposes.

It is well-equipped with brass or aluminum metal circular box, having a diameter of 100/125 mm.

A extremely magnetized needle on a hard steel pointed pivot is balanced at the centre of the box.

An Aluminum circle consists of a needle, which is graduated to 30 min. (0. 50) with engraved inverted figures.

The graduations are magnified by sliding the prism integrated with colored glasses having a sighting slit at the top.

The Prismatic Compass can be customized in conformance with the customers needs.

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)