

Product Name :
Prismatic Compass

Product Code :
NLAB-ENGINEERINGLB90001



Description :

Prismatic Compass

Technical Specification :

This precision engineered Prismatic Compass is one of the important devices which imparts immensely accurate readings for navigational purposes.

It is well-equipped with brass or aluminum metal circular box, having a diameter of 100/125 mm.

A extremely magnetized needle on a hard steel pointed pivot is balanced at the centre of the box.

An Aluminum circle consists of a needle, which is graduated to 30 min. (0. 50) with engraved inverted figures.

The graduations are magnified by sliding the prism integrated with colored glasses having a sighting slit at the top.

The Prismatic Compass can be customized in conformance with the customers needs.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com