


**Product Name :**  
Process Control Technology Lab Equipment

**Product Code :**  
NLAB-TECHNICALAB33001



**Description :**

Process Control Technology Lab Equipment

**Technical Specification :**

Capacitance Type

Digital Current Meter: Digital Indication of Current Output

Level Tank: Acrylic/Glass, PVC

Piping: Size ¼”

An english instruction manual consisting of experimental procedures, block diagram etc. will be provided along with the Apparatus.

The whole set-up is well designed and arranged on a rigid structure painted with industrial PU paint.

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com