


Product Name :
Proctor / CBR Test Kit According To NF Standard

Product Code :
NLAB-EDUCATIONTVET13023



Description :

Proctor / CBR Test Kit According To NF Standard

Technical Specification :

Normal Proctor mold according to NF

Modified Proctor / CBR mold complete with riser and perforated base, Ø152 x h152 mm

Dame Proctor normal according to NF

Dame Proctor modified according to NF

Leveling rule

CBR full base

Metallic fabric filter screen

Spacing disc with handle, h 25.4 mm

Spacing disc with handle, h 36 mm

Adjustable swelling plate

Comparator tripod

Comparator 25 x 0.01 mm

Ring overload 2300 g

Slotted overload 2300 g

Cutting kit

Filter paper, Ø150mm (100 pcs)

Immersion tank 600x400x400 mm

Universal extractor

Trowel

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003