


Product Name :
Proctor Compaction Test Apparatus

Product Code :
NLAB-TESTINGLABM109016



Description :

Proctor Compaction Test Apparatus

Technical Specification :

Soil compaction is of utmost requirement for constructing earth fill for dams, canal embankments, reservoirs, highways, runways and railways. The relationship between soil moisture content and compacted dry density largely determine the construction specifications and quality control of compacted earthfill. To check the strength of the soil the density of the soil is taken into account and hence the quality control is checked otherwise.

**Compaction Mould, complete with
Collar and Base Plate**

100 Proctor Compaction (27.3cm)
Height, 1000ml volume (Mild Steel)

225 Proctor Compaction (27.3cm)
Height, 2250ml volume (Mild Steel)

~~4381kg (Gun Metal)~~
~~4381kg (Gun Metal)~~
Height,2250ml volume (Gun Metal)

~~4381kg (Gun Metal)~~
~~4381kg (Gun Metal)~~
Height,1000ml volume (Gun Metal)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003