

**Product Name :**  
Proctor Compaction Test Apparatus

**Product Code :**  
NLAB-TESTINGLABM109016



**Description :**

Proctor Compaction Test Apparatus

**Technical Specification :**

Soil compaction is of utmost requirement for constructing earth fill for dams, canal embankments, reservoirs, highways, runways and railways. The relationship between soil moisture content and compacted dry density largely determine the construction specifications and quality control of compacted earthfill. To check the strength of the soil the density of the soil is taken into account and hence the quality control is checked otherwise.

**Compaction Mould, complete with Collar and Base Plate**

**100 Proctor Compaction (17.31) m**  
Height, 1000ml volume (Mild Steel)

**250 Proctor Compaction (17.31) m**  
Height, 2250ml volume (Mild Steel)

~~439~~ Brighton 400ml Dialer (27.8mm)  
Height, 2250ml volume (Gun Metal)

~~439~~ Brighton 400ml Dialer (27.8mm)  
Height, 1000ml volume (Gun Metal)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003