

Product Name :
Proctor Moulds and Rammers

Product Code :
NLAB-SOILELAB140032



Description :

Proctor Moulds and Rammers

Technical Specification :

Used for determining the relationship between the moisture content and density of compacted soil. The mould includes collar, mould body and base plate. The rammer construction includes a guide sleeve with vent holes.

Different versions are available that conform to the various commonly used standards. They are all made of plated steel and are identical in shape, only differing slightly in diameter and capacity.

In this section moulds and rammers conforming to ASTM/AASHTO are shown.

For the extrusion of soil specimens from the mould, the Universal specimen extruder may be used, see [Link to other products below](#).

An alternative (and preferable) method of compacting is to use an automatic compactor. For more information, see [Link to other products below](#).

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com