

**Product Name :**  
Proctor penetrometer

**Product Code :**  
NLAB-SOILELAB140021



**Description :**

Proctor penetrometer

**Technical Specification :**

Used for establishing the moisture content-penetration resistance relationship of fine-grained soils. It consists of a special spring dynamometer with pressure indicating scale on the stem of the handle. A sliding ring on the stem indicate the maximum pressure obtained in the test. Supplied in a wooden carrying case.

Load scale: 0 to 55 kg , 1 kg subdivisions with max load indicator

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003