

**Product Name :**  
Profoscope Type Pachometer

**Product Code :**  
NLAB-EDUCATIONTVET15032



**Description :**

Profoscope Type Pachometer

**Technical Specification :**

Location of metallic reinforcements and other metallic building elements and their orientation

Measurement of concrete cover thickness and reinforcement diameter

Functions for correcting the influence of adjacent bars

Compact and ergonomic indication device with display backlight

Assistance with optical and acoustic localization

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003