


Product Name : Programmable Dental Chair	Product Code : NLABMEDICAL-M32100030
--	--



Description :

Programmable Dental Chair

Technical Specification :

Dental Chair consists of a heavy duty programmable dental chair with two fixed and three programmable positions with synchronized motion of seat with backrest.

Standard Accessories of Dental Chair include:

- Two 3 Way Syringes
 - Purification system and integral spittoon
 - X-Ray Viewer
 - 2 points for high speed Air Turbine
 - One point for Air Motor
 - Automatic Cup Filler
 - Pneumatic Dental Stool
 - Water Spittoon and Saliva Ejector
 - Operating light mounted on smooth arms
- In addition to the standard accessories the inbuilt features of this dental unit include:
- Six point LED touch free sensor lighting unit
 - Rotatable Spittoon
 - LED illuminated Air Rotor Handpiece
 - Additional provision for operation of spittoon and tumbler from foot control

Inbuilt Peizoelectric Scaler with 5 tips

Micromotor with speed upto 36,000 RPM

Can be provided with ceramic white or coloured spittoon

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003