

Product Name : Programmable Dental Chair	Product Code : NLABMEDICAL-M32100030
Description : Programmable Dental Chair	
Technical Specification : Dental Chair consists of a heavy duty programmable dental chair with two fixed and three programmable positions with synchronized motion of seat with backrest. Standard Accessories of Dental Chair include: Two 3 Way Syringes Purification system and integral spittoon X-Ray Viewer 2 points for high speed Air Turbine One point for Air Motor Automatic Cup Filler Pneumatic Dental Stool Water Spittoon and Saliva Ejector Operating light mounted on smooth arms In addition to the standard accessories the inbuilt features of this dental unit include: Six point LED touch free sensor lighting unit Rotatable Spittoon LED illuminated Air Rotor Handpiece Additional provision for operation of spittoon and tumbler from foot control	

Inbuilt Piezoelectric Scaler with 5 tips

Micromotor with speed upto 36,000 RPM

Can be provided with ceramic white or coloured spittoon

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003