

Product Name :
Programmable Logic Controller 12 In 8 Out

Product Code :
NLAB-ENGINEERINGLB18002



Description :

Programmable Logic Controller 12 In 8 Out

Technical Specification :

Programmable Logic Controller 12 In 8 Out Lab Equipments Manufacturer, Suppliers & Exporters

Programmable Logic Controller “ 12 In/8 Out

Easy to be programmed from its own panel, without computer. It is composed (including the extension module) of 12 digital inputs and 8 relay outputs. Simulation of the inputs through switches and externally through terminals. The programming software is also supplied to allow programming the PLC from the computer, if so preferred

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

