

Product Name :
Programmable Logic Controller (MITSUBISHI PLC) Trainer

Product Code :
NLAB-ENGINEERINGLB45003



Description :

Programmable Logic Controller (MITSUBISHI PLC) Trainer

Technical Specification :

Programmable Logic Controller (MITSUBISHI PLC) Trainer Lab Equipments Manufacturer, Suppliers & Exporters Since PLC (Programmable Logic Controller) was firstly introduced in 1970 , it has been widely applied to various industrial uses such as machine and process controls. Designed with the latest microprocessor and electronic circuitry, today's compact-size PLCs feature high level of reliability, performance, speed and networking. The use of PLCs in automated production lines improves system reliability, product quality, information sharing, efficiency and flexibility and thus reduces costs. PLC-310 is a self-contained trainer which consists of a MITSUBISHI PLC main unit, I/O devices for simulation and I/O devices that are commonly used. It provides students with a thorough understanding of the theories and applications of programmable logic controllers. This trainer enables students to learn step by step the fundamentals of PLC and more advanced controls used in industry.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

