

Product Name :
Programmable Logic Controller Trainer

Product Code :
NLAB-EDUSCHOOLEQ189001



Description :

Programmable Logic Controller Trainer

Technical Specification :

The systems in our trainers are compatible with a variety of accessories and peripherals, thus allowing them to be easily expanded into larger systems.

The Programmable Logic Controller Trainer uses PLC products which are powerful and flexible systems ideally suited to your specialized needs.

A student experiment manual is provided to guide the student through the experiments.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003