

Product Name :
Programmable Logic Controller

Product Code :
NLAB-ENGINEERINGLB18001



Description :

Programmable Logic Controller 36 In 28 Out

Technical Specification :

Programmable Logic Controller 36 In 28 Out Lab Equipments Manufacturer, Suppliers & Exporters

1. The Programmable Logic Controller “36 In/28 Out unit is a programmable controller that combines high performances and ease of use for those who are approaching for the first time the world of PLC.

2. The proposed configuration includes:

CPU with 14 digital inputs, 10 relay outputs and 2 analogue inputs

Module with 16 digital inputs and 16 relay outputs

Module with 4 analogue inputs and 2 analogue outputs

3. The unit is fitted in an accessible container while on the panel the input/output terminals are shown, suitably duplicated by means of connectors

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003