

Product Name :
Projection Microscope

Product Code :
NLAB-EDUSCHOOLEQ173070

Description :

Projection Microscope

Technical Specification :

The microscopes are compact in design and can easily contrast and distortion free magnified images on a fine ground glass screen.

The structure of microscopic specimens becomes visible to the naked eye by taking the projection on the screen, during observation and demonstration to a group of students.

Nosepiece - Triple revolving nosepiece with accurate centering and positive click stops

Eyepieces - Huygenian 10X and 15X

Body - Monocular Fixed Straight Tube, length 160mm.

Screen - 150mm diameter Glass Projection Screen is provided, can be rotated over 360° on Microscope tube.

Stage - Mechanical stage of 135x120mm. Low co-axial drive controls for X-Y movement and cross movement 50X75mm, fully graduated vernier scale.

Illumination - 20W halogen lamp, 220/110Volts. Lamp slots provided for natural cooling of the illuminator.

Objectives - Achromatic 4X - 10X- 40XSL.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003