

Product Name :
Projection Slide Storage Box

Product Code :
NLAB-SCIENTIFICEQ15007



Description :

Projection Slide Storage Box

Technical Specification :

Useful to store 2x2? projection slides. Frame & slide holding grooves made of Teakwood duly numbered.

All fittings like hinges, clips, index card frames are of brass Lid clad in laminated sheet.

Capacity : 100 – 50 – 25 slides

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003