

Product Name :
Projection Slide Storage Cabinets

Product Code :
NLAB-SCIENTIFICEQ15008



Description :

Projection Slide Storage Cabinets

Technical Specification :

Projection Slide Storage Cabinet is available in:- For storage of 2x2? projection slides. Cabinet made of duro-board

clad in laminated sheet. Sliding drawers with wooden grooves carry slides vertically. An extra drawer is provided to keep index card sheets. Dust proof door has locking arrangement.

Capacity 600 slides with 5 drawers of 120 slides each.

Capacity 1000 slides with 5 drawers of 200 slides each.

Capacity 2000 slides with 10 drawers of 200 slides each.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

