

Product Name :
Projections Epidiascope

Product Code :
NLAB-SCIENTIFICEQ15009



Description :

Projections
Epidiascope

Technical Specification :

Epidiascope is a convertible projection equipment with 'A' grade anastigmatic lenses and surface-coated reflecting mirrors for brilliant projections of opaque objects as well as projection slides, film strips etc.

Specifications:

- Noise-free cooling system for long time projections.
- Complete with attachments for projection of 2 x 2? projection slides, 35mm film strips and prepared micro-slides.
- Provided with 100W-220V ac halogen bulb, cord & plug.
- Fit for sharp projection upto 30 ft.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003