

**Product Name :**  
Projections Epidiascope

**Product Code :**  
NLAB-SCIENTIFICEQ15009



**Description :**

Projections  
Epidiascope

**Technical Specification :**

Epidiascope is a convertible projection equipment with 'A' grade anastigmatic lenses and surface-coated reflecting mirrors for brilliant projections of opaque objects as well as projection slides, film strips etc.

**Specifications:**

- Noise-free cooling system for long time projections.
- Complete with attachments for projection of 2 x 2? projection slides, 35mm film strips and prepared micro-slides.
- Provided with 100W-220V ac halogen bulb, cord & plug.
- Fit for sharp projection upto 30 ft.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003