

Product Name :
Pull Off Adhesion Tester

Product Code :
NLAB -CONCRLAB120045



Description :

Pull Off Adhesion Tester

Technical Specification :

Self-Alignment Feature

Unequal pulling force during testing caused by uneven adhesive bond lines and coating surfaces can result in random, unexplainable readings. To obtain more repeatable and meaningful adhesion measurements, it is imperative that the pulling force applied to the test dolly is uniformly distributed over the surface being tested. The Self-Aligning Dolly enables uniform distribution of the pulling force over the surface being tested, preventing a one-sided pull-off.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003