

**Product Name :**  
Pulley Mounted Rod

**Product Code :**  
NLAB-EDUSCHOOLEQ1740106



**Description :**

Pulley Mounted Rod

**Technical Specification :**

Single 50mm diameter ball bearing pulley, mounted in aluminium boss provided with a hole and clamping screw to accommodate rod.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003