

Product Name :
Pulse Oximeter System

Product Code :
NLABMEDICAL-M3410001



Description :

Pulse Oximeter System

Technical Specification :

Large screen: SpO₂,PR are shown on LED screen, Pleth waveform is displayed on 3.5" LCD screen in real time

Light in weight, only 1kg with battery

72-hour data storage, 100-item patient ID

72-hour SpO₂ and PR spot graphs

3-level visual & audio alarmIntellectual control of power consumption, internal rechargeable battery lasts 24 hours or AC power

Battery-low indication

Built-in thermal recorder(optional)

Suitable for adult, pediatric and neonatal

Widely used in small clinic, all departments in hospital, ambulance

Compatible with CHOICE central monitoring system

USB data upload (optional)

Anti-high-frequency electrosurgical equipment

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003